



Product Datasheet

DHRS13 antibody (orb354583)

Description	DHRS13 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Dehydrogenase/reductase SDR family member 13 protein (1-327AA)
Target	DHRS13
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6UX07
Dilution Range	WB:1:500-1:5000, IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293 whole ce...



Immunohistochemical staining of human br...