

Product Datasheet

DHX30 antibody (FITC) (orb354597)

Description	DHX30 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Putative ATP-dependent RNA helicase DHX30 protein (694-1000AA)
Target	DHX30
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q7L2E3
Expiration Date	12 months from date of receipt.