

Product Datasheet

DNAAF5 antibody (orb354623)

Description	DNAAF5 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Dynein assembly factor 5, axonemal protein (621-843AA)
Target	DNAAF5
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q86Y56
Expiration Date	12 months from date of receipt.