

Product Datasheet

DYNLL2 antibody (HRP) (orb354686)

Description	DYNLL2 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Dynein light chain 2, cytoplasmic protein (1-89AA)
Target	DYNLL2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96FJ2
Expiration Date	12 months from date of receipt.