

Product Datasheet

FBXW10 antibody (Biotin) (orb354916)

Description	FBXW10 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human F-box/WD repeat-containing protein 10 protein (701-1000AA)
Target	FBXW10
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5XX13
Expiration Date	12 months from date of receipt.