

Product Datasheet

FXYP7 antibody (Biotin) (orb354980)

Description	FXYP7 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human FXYP domain-containing ion transport regulator 7 protein (47-80AA)
Target	FXYP7
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P58549
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.