

Product Datasheet

CD300LB antibody (orb35500)

Description	Rabbit polyclonal antibody to CD300LB (CD300 molecule-like family member b), also known as CLM7, IRE...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	CD300LB
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	22689
Uniprot ID	A8K4G0
NCBI	NP_777552.3
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

