

Product Datasheet

SPRR2A antibody (orb35502)

Description	SPRR2A antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Small proline-rich protein 2A protein (1-72AA)
Target	SPRR2A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P35326
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.