


Product Datasheet

GNGT2 antibody (orb355079)

Description	GNGT2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Recombinant Human Guanine nucleotide-binding protein G(I)/G(S)/G(O) subunit gamma-T2 protein (1-66AA)
Target	GNGT2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O14610
Dilution Range	IF:1:200-500
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of A549 cell...