



Product Datasheet

GPR18 antibody (orb355103)


Description	GPR18 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human N-arachidonyl glycine receptor protein (290-331AA)
Target	GPR18
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q14330
Dilution Range	WB:1:500-2000, IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U251 cell...



Immunohistochemical staining of human te...



Western blot analysis of Mouse brain tis...