

Product Datasheet

CMTM5 antibody (FITC) (orb35512)

Description	CMTM5 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human CKLF-like MARVEL transmembrane domain-containing protein 5 protein (1-35AA)
Target	CMTM5
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96DZ9
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.