

Product Datasheet

HDHD3 antibody (Biotin) (orb355164)

Description	HDHD3 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Haloacid dehalogenase-like hydrolase domain-containing protein 3 protein (1-251AA)
Target	HDHD3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9BSH5
Expiration Date	12 months from date of receipt.