

Product Datasheet

HIF1AN antibody (HRP) (orb355186)

Description	HIF1AN antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Hypoxia-inducible factor 1-alpha inhibitor protein (2-254AA)
Target	HIF1AN
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NWT6
Expiration Date	12 months from date of receipt.