



Product Datasheet

HM13 antibody (orb355195)

Description	HM13 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant Human Minor histocompatibility antigen H13 protein (336-377AA)
Target	HM13
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8TCT9
Dilution Range	WB:1:500-2000
Expiration Date	12 months from date of receipt.



Immunofluorescent analysis of 293 cells ...



Western blot analysis of Mouse liver ti...