

Product Datasheet

HMGCLL1 antibody (Biotin) (orb355200)

Description	HMGCLL1 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human 3-hydroxymethyl-3-methylglutaryl-CoA lyase, cytoplasmic protein (2-340AA)
Target	HMGCLL1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8TB92
Expiration Date	12 months from date of receipt.