



Product Datasheet

HSD17B12 antibody (orb355235)


Description	HSD17B12 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Very-long-chain 3-oxoacyl-CoA reductase protein (40-181AA)
Target	HSD17B12
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q53GQ0
Dilution Range	WB:1:500-2000, IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of U251(Lane 1), A...



Immunohistochemical staining of human br...



Immunofluorescence analysis of U251 cell...