


Product Datasheet

Hspa1b antibody (orb355243)

Description	Hspa1b antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Mouse Heat shock 70 kDa protein 1B protein (2-642AA)
Target	Hspa1b
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P17879
Dilution Range	WB:1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of A431(Lane 1), A...