




Product Datasheet

IL5RA antibody (orb35532)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthetic Peptide
Target	IL5RA
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	47685
Uniprot ID	Q01344
NCBI	NP_783853.1 , NP_783852.1 , NP_783851.1 , NP_783855.1 , NP_000555.2 , NP_783854.1 , NP_001230028.1
Dilution Range	WB: 1:1000, IHC-P: 1:50-100, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of IL5RA antibody...



Immunohistochemical analysis of formalin-fixed tissue...