

Product Datasheet

KRT33B antibody (Biotin) (orb355396)

Description	KRT33B antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Keratin, type I cuticular Ha3-II protein (1-404AA)
Target	KRT33B
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q14525
Expiration Date	12 months from date of receipt.