

Product Datasheet

L3HYPDH antibody (Biotin) (orb355404)

Description	L3HYPDH antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Trans-3-hydroxy-L-proline dehydratase protein (13-312AA)
Target	L3HYPDH
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96EM0
Expiration Date	12 months from date of receipt.