

Product Datasheet

LECT1 antibody (Biotin) (orb355420)

Description	LECT1 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Leukocyte cell-derived chemotaxin 1 protein (215-334AA)
Target	CNMD
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	075829
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.