



Product Datasheet

LGALS12 antibody (orb355423)


Description	LGALS12 antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Galectin-12 protein (1-336AA)
Target	LGALS12
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96DT0
Dilution Range	WB:1:500-2000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human brain tissue.



Immunohistochemical staining of human skin tissue.



Western blot analysis of A549 whole cell lysate and kidney tissue.