

Product Datasheet

Ly6e antibody (FITC) (orb355501)

Description	Ly6e antibody (FITC)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	FITC
Immunogen	Recombinant Mouse Lymphocyte antigen 6E protein (21-102AA)
Target	Ly6e
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q64253
Expiration Date	12 months from date of receipt.