



Product Datasheet

MAS1L antibody (orb355571)


Description	MAS1L antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Mas-related G-protein coupled receptor MRG protein (318-378AA)
Target	MAS1L
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P35410
Dilution Range	WB:1:500-2000, IHC-P:1:20-200, IF:1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of colo320(Lane 1)...



Immunofluorescence analysis of HepG2 cel...



Immunofluorescence analysis of HepG2 cel...