

Product Datasheet

MIEF1 antibody (orb355619)

Description	MIEF1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Mitochondrial dynamics protein MID51 protein (51-309AA)
Target	MIEF1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NQG6
Dilution Range	WB:1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse liver tis...