

Product Datasheet

MIF4GD antibody (Biotin) (orb355628)

Description	MIF4GD antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human MIF4G domain-containing protein (1-133AA)
Target	MIF4GD
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	A9UHW6
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.