

Product Datasheet

MRPL18 antibody (HRP) (orb355658)

Description	MRPL18 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human 39S ribosomal protein L18, mitochondrial protein (14-180AA)
Target	MRPL18
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9H0U6
Expiration Date	12 months from date of receipt.