

Product Datasheet

GLT1D1 antibody (FITC) (orb35567)

Description	GLT1D1 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Glycosyltransferase 1 domain-containing protein 1 protein (17-152AA)
Target	GLT1D1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96MS3
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.