

Product Datasheet

SPG16 antibody (orb35571)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Synthetic Peptide
Target	SPAG16
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	70818
Uniprot ID	Q8N0X2
NCBI	NP_001020607.1, NP_078808.3
Dilution Range	IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.

