

Product Datasheet

MRRF antibody (Biotin) (orb355728)

Description	MRRF antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Ribosome-recycling factor, mitochondrial protein (56-262AA)
Target	MRRF
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96E11
Expiration Date	12 months from date of receipt.