



Product Datasheet

MTR antibody (orb355763)


Description	MTR antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP
Immunogen	Recombinant Human Methionine synthase protein (923-1265AA)
Target	MTR
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q99707
Dilution Range	IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human adrenal gland tissue.



Immunoprecipitation analysis of HeLa cell lysates.



Immunohistochemical staining of human placental tissue.