

Product Datasheet

MYNN antibody (FITC) (orb355781)

Description	MYNN antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Myoneurin protein (51-330AA)
Target	MYNN
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NPC7
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.