

Product Datasheet

NKIRAS1 antibody (FITC) (orb355837)

Description	NKIRAS1 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human NF-kappa-B inhibitor-interacting Ras-like protein 1 protein (1-192AA)
Target	NKIRAS1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NYS0
Expiration Date	12 months from date of receipt.