


Product Datasheet

NUDCD3 antibody (orb355895)

Description	NUDCD3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human NudC domain-containing protein 3 protein (142-361AA)
Target	NUDCD3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q81VD9
Dilution Range	WB:1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of 293T(Lane 1), H...