

Product Datasheet

ORMDL1 antibody (Biotin) (orb355932)

Description	ORMDL1 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human ORM1-like protein 1 protein (48-100AA)
Target	ORMDL1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9P0S3
Expiration Date	12 months from date of receipt.