

Product Datasheet

PDIK1L antibody (FITC) (orb356025)

Description	PDIK1L antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Serine/threonine-protein kinase PDIK1L protein (1-341AA)
Target	PDIK1L
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8N165
Expiration Date	12 months from date of receipt.