

Product Datasheet

NDUFB5 antibody (HRP) (orb356046)

Description	NDUFB5 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Prefoldin subunit 4 protein (2-133AA)
Target	PFDN4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NQP4
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.