

Product Datasheet

PIH1D1 antibody (Biotin) (orb356072)

Description	PIH1D1 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human PIH1 domain-containing protein 1 protein (1-290AA)
Target	PIH1D1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NWS0
Expiration Date	12 months from date of receipt.