

Product Datasheet

Plet1 antibody (Biotin) (orb356112)

Description	Plet1 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Mouse Placenta-expressed transcript 1 protein (28-218AA)
Target	Plet1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8VEN2
Expiration Date	12 months from date of receipt.