


Product Datasheet

POLR3F antibody (orb356147)

Description	POLR3F antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human DNA-directed RNA polymerase III subunit RPC6 protein (1-270AA)
Target	POLR3F
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9H1D9
Dilution Range	WB:1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa (Lane 1), H...