



Product Datasheet

PPP2R2B antibody (orb356187)


Description	PPP2R2B antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Immunogen	Recombinant Human Serine/threonine-protein phosphatase 2A 55 kDa regulatory subunit B beta isoform protein (1-200AA)
Target	PPP2R2B
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q00005
Dilution Range	IHC-P:1:50-200, IF:1:50-500
Expiration Date	12 months from date of receipt.



Western blot analysis of SH-SY5Y(Lane 1)...



Immunohistochemical staining of human CO...



Immunofluorescence analysis of SH-SY5Y C...