



Product Datasheet

PPP2R4 antibody (orb356191)


Description	PPP2R4 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Serine/threonine-protein phosphatase 2A activator protein (1-250AA)
Target	PTPA
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q15257
Dilution Range	WB:1:500-5000, IHC-P:1:1000-2000, IF:1:50-500
Expiration Date	12 months from date of receipt.



Western blot analysis of K562(Lane 1), M...



Immunohistochemical staining of human pr...



Immunohistochemical staining of human en...