

Product Datasheet

PRKRA antibody (Biotin) (orb356212)

Description	PRKRA antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Interferon-inducible double-stranded RNA-dependent protein kinase activator A protein (1-313AA)
Target	PRKRA
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	075569
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.