

Product Datasheet

Psat1 antibody (FITC) (orb356241)

Description	Psat1 antibody (FITC)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	FITC
Immunogen	Recombinant Mouse Phosphoserine aminotransferase protein (1-370AA)
Target	Psat1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q99K85
Expiration Date	12 months from date of receipt.