



Product Datasheet

RAB11B antibody (orb356279)

Description	RAB11B antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Ras-related protein Rab-11B protein (2-215AA)
Target	RAB11B
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q15907
Dilution Range	WB:1:500-2000, IHC-P:1:20-200
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa cells (Lane 1 and Lane 2)



Immunohistochemical staining of human placenta (pl...)