

Product Datasheet

RAB37 antibody (FITC) (orb356293)

Description	RAB37 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Ras-related protein Rab-37 protein (2-220AA)
Target	RAB37
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96AX2
Expiration Date	12 months from date of receipt.