


Product Datasheet

RBP5 antibody (orb356363)

Description	RBP5 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Retinol-binding protein 5 protein (1-135AA)
Target	RBP5
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P82980
Dilution Range	WB:1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse liver tis...