



Product Datasheet

RPL10L antibody (orb356483)


Description	RPL10L antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human 60S ribosomal protein L10-like protein (1-214AA)
Target	RPL10L
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96L21
Dilution Range	WB:1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat(Lane 1),...



Immunohistochemical staining of human me...



Immunohistochemical staining of human ce...