

Product Datasheet

SEC31A antibody (FITC) (orb356645)

Description	SEC31A antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Protein transport protein Sec31A protein (335-503AA)
Target	SEC31A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O94979
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.