


Product Datasheet

SIGLEC15 antibody (orb356703)

Description	SIGLEC15 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Sialic acid-binding Ig-like lectin 15 protein (20-263AA)
Target	SIGLEC15
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6ZMC9
Dilution Range	WB:1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse spleen ti...